

ProductProfile

Hematology Controls

Avantor™ offers different types of J.T.Baker® hematology controls, suitable for use with many popular hematology instruments. Below is an overview for what instrument the controls can be used.

J.T.Baker® Hematology Controls Selection Guide

Instrument	Type	8 PMC	8 PMC VET	3-Diff	5D	ADV-Diff	BC-Diff 5	CD-Diff	K-Diff	SF-Diff	XE-Diff	Platelet	Cal-Set 1
Abbott®	Cell-Dyn® Sapphire, Ruby, 4000, 3700, 3500							x				x	
	Cell Dyn® 3200 and 3000							x					
	Cell-Dyn® 1800, 1700, 1600 and 1300	x		x									x
	Cell-Dyn® Emerald 18			x									
ABX®	Pentra® 80, 60, 60C, 60C+, 60ES, 60MS	x										x	x
	Micros® 60, ES60, 45	x		x								x	x
Avantor Performance Materials™	BeneSphera™ H34, H32 3-part diff	x		x									x
	BeneSphera™ H32 3-part diff VET		x										
AVL	818, 816	x		x									x
Bayer®	Advia® 60	x		x									x
	Advia® 120					x						x	
Beckman Coulter®	AcT Diff™, AcT Diff 2™	x		x									x
	AcT10™, AcT8™, AcT5™	x											x
	STKS™, GenS™, LH780, LH750, LH500, HmX™, MaxM™	x			x								
Diatron	Abacus Junior 5									x			x
	Abacus, Abacus Junior	x		x									x
Drew	D3, Excell 18	x		x									
Erma	PCE-210, PCE-170	x		x									x
HTI®	Micros CC20, CC18			x									
Medonic	CA620-20, CA620-16, CA610, CA600, CA570, CA530-16			x									x
	CA590, CA580, CA530-9	x											x
Mindray®	BC-5380, BC-5300						x						x
	BC-3200, BC-3000+, BC-2800, BC-2300, BC-2000	x		x									x
Nihon Kohden®	Celltac™ F MEK-8222K, α MEK-8118K+QA-810V, ES MEK-7300, E MEK-7222K	x						x					x
	Celltac™ α MEK-6450, MEK-6420, MEK-6410, MEK-6400, MEK-6318K, MEK-6108	x		x									x
Orphee	Mythic™ 18			x									x
Rayto	RT-7600			x									x
SEAC	H20 Genius	x						x					
	H12, H8	x											x
Swelab®	AC970eo+, AC920eo+	x										x	
	AC910, AC900											x	x
Sysmex®	SF3000									x			x
	K4500, K1000, KX21, KX21-N	x							x				x
	XE2100, XT2000i, XT1800i, XS1000i										x	x	
Urit®	3010, 2900			x									

The information contained herein has not been approved by listed instrument manufacturer. Always refer to the user manual provided with the equipment at issue.

Stability

To guarantee the stability in the table below, the hematology controls have to be stored upright at 2-8 °C when not in use. The vials may be stressed for 48 hours at 18 °C maximum, without losing product integrity. The stability period of opened vials is guaranteed when the hematology controls are used in accordance with the instructions for use.

J.T.Baker® Hematology Controls Overview

J.T.Baker® Product	Level	Control Number	Volume	Pack type	Stability (months)	Open vial stability (weeks)	Ordering frequency (per year)
8-Parameter Control	L/N/H	3427/3428/3429	2.5 ml	Screw cap	4	4	4
	L/N/H	3463/3464/3465	2.5 ml	Pierceable cap	4	4	4
	L/N/H	3701/3702/3703	4.5 ml	Pierceable cap	4	4	4
	L+N+H	3746	3 x 2.5 ml	Screw cap	4	4	4
	4xN	3747	4 x 2.5 ml	Screw cap	4	4	4
	1xL+4xN+1xH	3751	6 x 2.5 ml	Screw cap	4	4	4
8-Parameter Control extended	L/N/H	3633/3634/3635	2.5 ml	Screw cap	6	4	4
8-Parameter Control VET	L/N/H	2904/2905/2906	2.5 ml	Screw cap	4	4	4
3-Diff Control	L/N/H	3433/3434/3435	2.5 ml	Screw cap	4	3	4
	L/N/H	3502/3503/3504	4.5 ml	Pierceable cap	4	3	4
	4xL/N/H	3466/3467/3468	4 x 2.5 ml	Screw cap	4	3	4
3-Diff Control extended	L/N/H	3421/3422/3423	2.5 ml	Screw cap	6	3	4
5D Control	L/N/H	3681/3682/3683	5.0 ml	Pierceable cap	3	2	6
ADV-Diff Control	L/N/H	3684/3685/3686	3.5 ml	Screw cap	3	2	6
BC-Diff 5 Control	L/N/H	3613/3614/3615	3.0 ml	Pierceable cap	3	2	6
CD-Diff Control	L/N/H	3452/3453/3454	3.0 ml	Pierceable cap	3	2	6
	2xL+2xN+2xH	3838	6 x 3.0 ml	Pierceable cap	3	2	6
K-Diff Control	L/N/H	3455/3456/3457	2.5 ml	Screw cap	4	3	4
SF-Diff Control	L/N/H	3693/3694/3695	4.5 ml	Pierceable cap	3	2	6
XE-Diff Control	L/N/H	3731/3732/3733	4.5 ml	Pierceable cap	3	2	6
Cal-Set 1		3940	2 x 2.5 ml	Screw cap	3	2	4
Platelet - Extended value		3424	5 x 3 ml	Pierceable cap	3	2	6
WBC Reduced RBC	L/H	3698/3699	3.0 ml	Pierceable cap	3	4	6

Abbott® and Cell-Dyn® are trademarks of Abbott Laboratories Corporation. ABX®, Pentra® and Micros® are trademarks of HORIBA ABX SAS. AVL is manufactured by Roche Diagnostics Corporation. Bayer® and Advia® are trademarks of Bayer Corporation. Beckman Coulter®, ACT Diff®, ACT Diff 2®, ACT10®, ACT8®, ACT5®, STKS®, GenS®, LH780, LH750, LH500, HmX®, MaxM® are trademarks of Beckman Coulter, Inc. Diatron is manufactured by Diatron MI PLC. Drew is manufactured by Drew Scientific, Ltd. Erma is manufactured by Erma, Inc. Hospitex is manufactured by Hospitex Diagnostics. HTT® is a trademark of High Technology, Inc. Medonic is manufactured by Boule Medical AB. S.A. Mindray® is a trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. Nihon Kohden® is a trademark of Nihon Kohden Corporation. Orphee S.A. is manufactured by The Cormay Group, Mythic™ is trademark of Orphee S.A. Rayto is manufactured by Rayto Life and Analytical Sciences Co., Ltd. Seac is manufactured by RADIM Spa. Swelab® is a trademark of Boule Medical AB. Sysmex® is a trademark of Sysmex Corporation. Urit® is manufactured by Urit Medical Electronic Co., Ltd.



Phillipsburg, NJ 9001:2008 & 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, the Netherlands 9001:2008, 14001:2004 & 13485:2003

Gliwice, Poland 9001:2008 & 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India 9001:2008, 14001:2004 & 13485:2003
Panoli, India 9001:2008



Ordering Information and Assistance

Customer Service and Technical Service

India

1201-1206, Pinnacle Business Tower
Shooting Range Road, Surajkund, Faridabad-121009, India
TEL: +91-129-4267000, 4267100, 4267200
FAX: +91-129-4267299, 4267199
E-MAIL: diagnostics.india@avantormaterials.com

Netherlands

Teugseweg 20, 7418 AM Deventer, the Netherlands
TEL: +31-570-687500
FAX: +31-570-687574
E-MAIL: avantor.emea@avantormaterials.com
www.avantormaterials.com

ASK Avantor™

Our Web site features ASK Avantor™, which includes live chat capabilities with customer service representatives.
www.avantormaterials.com/askavantor

Worldwide Locations

• China • Malaysia • North America
• India • Mexico • Poland
• Korea • Netherlands • Taiwan

For contact information at these locations, visit
www.avantormaterials.com/WorldwideDirectory

