


BeneSphera™

# 3 PART DIFFERENTIAL

Hematology Analyzer



## Features and Benefits

- **3 Part differentiation of WBC** – 20 parameters + 3 histograms
- **2 Counting modes** – whole blood and prediluted
- **Throughput** – 60 samples per hour
- **Touch screen** – with 10.4" TFT color LCD
- **Presentation of results** – built-in thermal recorder, optional external printer
- **Sample results** – storage up to 50,000 sample results (including histograms)
- **Multi-language software** – available on request
- **Dedicated J.T.Baker® reagents** – providing reliable and consistent results
- **Regulatory compliance with certifications** – ISO and CE Certified
- **Quality control** – to confirm the precision and accuracy of instrument operation
- **Worldwide presence** – global manufacturing capability and sales network



BeneSphera™



Part of



Sample Test

Auto Calibration Setting



Manual Calibration



L-J QC Setting



Profile

## Technical Specifications

<b>Product</b>	3 Part Differential Hematology Analyzer	
<b>Counting Channel</b>	2	
<b>Throughput</b>	60 tests/hour	
<b>Parameters</b>	20	
<b>Operating</b>	Touch Screen, Mouse, Keyboard	
<b>Sample Volume</b>	Whole blood 9.8 µl, Anti-coagulated Peripheral 9.8 µl, Prediluted 20µl	
<b>Reagents</b>	Diluent (15.5ml/test), Lyse (0.6ml/test), Cleanser (1.2 ml/test)	
<b>Carry Over</b>	WBC, RBC, HGB < 0.5%, PLT < 1%	
<b>QC Mode</b>	XB & L-J; X-R; Westgard Rules	
<b>Storage</b>	50,000	
<b>Interface</b>	RS232C, USB*2, Keyboard, Mouse Interface (PS2), LAN	
<b>Linearity Range</b>	WBC (10 <sup>9</sup> /L) 0 - 99.9 RBC (10 <sup>12</sup> /L) 0 - 9.99 HGB (g/L) 0 - 300 MCV (fL) 40 - 150 PLT (10 <sup>9</sup> /L) 0 - 999	
<b>Precision (CV%)</b>	WBC (10 <sup>9</sup> /L) ≤ 2.0% (7.0 - 15*10 <sup>9</sup> /L) RBC (10 <sup>12</sup> /L) ≤ 1.5% (3.5 - 6*10 <sup>12</sup> /L) HGB (g/L) ≤ 1.5% (110 - 180 g/L) MCV (fL) ≤ 0.4% (80 - 110 FL) PLT (10 <sup>9</sup> /L) ≤ 4.0% (100 - 500*10 <sup>9</sup> /L)	
<b>Power Requirement</b>	AC 110V~220V±10%, 50-60Hz, External AC Adaptor	
<b>Dimension</b>	436 x 363 x 367 mm	
<b>Weight</b>	18 Kg	



Phillipsburg, NJ 9001:2008 & 14001:2004  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, the Netherlands 9001:2008, 14001:2004 & 13485:2003

Gliwice, Poland 9001:2008 & 17025:2005  
Selangor, Malaysia 9001:2008  
Dehradun, India 9001:2008, 14001:2004 & 13485:2003  
Mumbai, India 9001:2008 & 17025:2005



## About Avantor™ Performance Materials

Avantor Performance Materials manufactures and markets high-performance chemistries and materials around the world under several respected brand names, including the J.T.Baker®, Macron Fine Chemicals™, Rankem™, BeneSphera™ and POCH™ brands.

Avantor products are used in a wide range of industries. Our biomedical and life science solutions are used in pharmaceutical production, laboratory research for academic, industry and quality control, and in medical lab testing. Our electronics materials products are used in the manufacturing of semiconductors. Based in Center Valley, Pennsylvania (USA), Avantor is owned by an affiliate of New Mountain Capital, LLC.

For additional information please visit [www.avantormaterials.com](http://www.avantormaterials.com) or follow [www.twitter.com/avantor\\_news](https://twitter.com/avantor_news)

## Ordering Information and Assistance

### Customer Service and Technical Service India

1201-1206, Pinnacle Business Tower  
Shooting Range Road, Surajkund, Faridabad - 121009, India  
TEL: +91-129-4267000, 4267100, 4267200  
FAX: +91-129-4267299, 4267199  
E-MAIL: [diagnostics.india@avantormaterials.com](mailto:diagnostics.india@avantormaterials.com)

### Netherlands

Teugseweg 20, 7418 AM Deventer, the Netherlands  
TEL: +31-570-687500  
FAX: +31-570-687574  
E-MAIL: [avantor.emea@avantormaterials.com](mailto:avantor.emea@avantormaterials.com)

[www.avantormaterials.com](http://www.avantormaterials.com)

### ASK Avantor™

Our Web site features ASK Avantor™ which includes live chat capabilities with customer service representatives.  
[www.avantormaterials.com/askavantor](http://www.avantormaterials.com/askavantor)

## Corporate Headquarters

Avantor Performance Materials, Inc.  
3477 Corporate Parkway  
Suite #200  
Center Valley, PA 18034  
USA

## Worldwide Locations

- China
- India
- Korea
- Malaysia
- Mexico
- Netherlands
- North America
- Poland
- Taiwan

For contact information at these locations, visit [www.avantormaterials.com/WorldwideDirectory](http://www.avantormaterials.com/WorldwideDirectory)