

# iM50

## Patient Monitor

Version 1.4



## Main Unit Specification

### Physical Specifications

Dimension	261 mm (W) × 215 mm (H) × 198 mm (D)
Weight	< 3.6 kg (standard configuration, without battery)

### Power Supply

Power Supply	100 V to 240 V~, 50 Hz/60 Hz
Current	1.0 A-0.5 A

### Battery

Battery Type	Rechargeable lithium-ion battery	
Capacitance	2500 mAh, 5000 mAh	
Operating Time	2500 mAh	≥4.5 h
	5000 mAh	≥10 h
Charge Time	2500 mAh	≤3.5 h, 100% charge
		≤3.15 h, 90% charge
	5000 mAh	≤6.5 h, 100% charge
		≤5.85 h, 90% charge

### Display

Display screen	8.4 inch color TFT, supporting touch screen
Resolution	800 × 600
Wave	A maximum of 13 waveforms

### Recorder

Record Width	48 mm
Paper Speed	12.5 mm/s, 25 mm/s, 50 mm/s
Channels	3
Recording types	Continual real-time recording
	8-second real-time recording
	20-second real-time recording,
	Trend graph recording
	Trend table recording
	NIBP review recording
	Arrhythmia review recording
	Alarm review recording
	Drug calculation titration recording
	Hemodynamic Calculation result recording
	Oxygenation Calculation result recording
	Ventilation Calculation result recording
	Renal Function Calculation result recording
	12-lead diagnosis recording
	Frozen waveform recording
	ST view recording
	QT view recording

### Data Storage

#### Internal Temporary Memory

A single piece of patient data maximally contains the following information:

Trend Data	3 hour, at 1 s resolution
	120 hours, at 1 min resolution

Alarm Events	Up to 200 sets
NIBP Measurement Data	1200 sets
Arrhythmia Events	Up to 200 sets
12-lead Analysis Result	Up to 50 sets

Full disclosure waveforms	48 hours
---------------------------	----------

#### Storage capacity for 1G extended space:

Continuous parameter data	720 hours, resolution: 1 min
NIBP data	At least 68000 sets
Physiological alarm event	At least 4500 sets
Arrhythmia event	At least 4500 sets
Full disclosure waveforms	30 hours

### Wi-Fi

IEEE	802.11a/b/g/n
Frequency Band	2.4 GHz ISM band & 5 G ISM band

### Interfaces and others

Nurse Call / Analog Output/ Defibrillator Synchronization	1
USB Interfaces	2
VGA Interface	1
Network Interface	1
Anti-theft lock interface	1

### ECG

Lead Mode	3 Electrodes: I, II, III
	5 Electrodes: I, II, III, aVR, aVL, aVF, V
	6 Electrodes: I, II, III, aVR, aVL, aVF, Va, Vb.
	10 Electrodes: I, II, III, aVR, aVL, aVF, V1-V6
Electrode Standard	AHA, IEC
Gain Selection	×0.125, ×0.25, ×0.5, ×1, ×2, ×4, AUTO gain
Sweep	6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s
Bandwidth (-3 dB)	Diagnosis: 0.05 Hz to 150 Hz
	Diagnosis 1: 0.05 Hz to 40 Hz
	Monitor: 0.5 Hz to 40 Hz
	Surgery: 1 Hz to 20 Hz
	Enhanced: 2 Hz ~18 Hz
CMRR	Customized: High-pass Filter and Low-pass Filter
	Diagnosis: > 95 dB
	Diagnosis 1: > 105 dB (when Notch is turned on)

	Monitor: > 105 dB Surgery: > 105 dB Enhanced: > 105 dB Customized: > 105 dB (Low-pass Filter < 40 Hz) > 95 dB (Low-pass Filter > 40 Hz)
<b>Hum Filter</b>	In diagnosis, Diagnosis 1, monitor, surgery, enhanced and customized modes: 50 Hz/60 Hz (Hum Filter can be turned on or off manually)
<b>Recovery Time After Defibrillation</b>	< 5 s (measured without electrodes as IEC60601-2-27:2011, Sect. 201.8.5.5.1 requires.)
<b>ESU Protection</b>	Cut mode: 300 W Coagulation mode: 100 W Restore time: ≤10 s
<b>Pace Pulse Detecting Lead</b>	one among I, II, III, aVR, aVL, aVF, V1-V6

<b>Heart Rate</b>	
<b>Range</b>	ADU: 15 bpm to 300 bpm PED/NEO: 15 bpm to 350 bpm
<b>Accuracy</b>	±1% or ±1 bpm, whichever is greater
<b>Resolution</b>	1 bpm

<b>PVC</b>	
<b>Range</b>	ADU: (0 to 300) PVCs/ min PED/NEO: (0 to 350) PVCs/ min
<b>Resolution</b>	1 PVCs/min

<b>ST value</b>	
<b>Range</b>	-2.0 mV to +2.0 mV
<b>Accuracy</b>	-0.8 mV to +0.8 mV: ±0.02 mV or 10%, whichever is greater. Beyond this range: not specified.
<b>Resolution</b>	0.01 mV

<b>QT measurement</b>	
<b>Range</b>	200 ms ~ 800 ms
<b>Resolution</b>	4 ms
<b>Accuracy</b>	± 30 ms

<b>QTc measurement</b>	
<b>Range</b>	200ms ~ 800 ms
<b>Resolution</b>	1 ms

<b>ΔQTc measurement</b>	
<b>Range</b>	-600 ms ~ 600 ms
<b>Resolution</b>	1 ms

<b>Arrhythmia analysis</b>	
Asystole, Sustain VT, V-Fib/V-Tach, ExtremeTachy, ExtremeBrady, V-Tach, Vent Brady, Tachy, Brady, Wide QRS Tachy, Non-Sustain VT, Afib, Vent Rhythm, Acc. Vent Rhythm, Pause, Pauses/min High, PVCs High, R on T, PVC Bigeminy, PVC Trigeminy, Pacer not Pacing, Pacer not Capture, Missed Beat, VEB, PVC, Couplet, Run PVCs, IPVC, Irr Rhythm, PAC Bigeminy, Multiform PVCs, PAC Trigeminy, Low Voltage (Limb)	

<b>12-Lead ECG Synchronization Analysis</b>	
Average parameters of heart beat	PR interval (ms)
Heart rate (bpm)	QRS interval (ms)
Time limit of P wave (ms)	QT/QTc (ms)
P-QRS-T AXIS	

<b>RESP</b>	
<b>Method</b>	Impedance between RA-LL, RA-LA

<b>Measurement lead</b>	Options are lead I and II. The default is lead II.
<b>RR Measuring Range</b>	Adult: 0 rpm to 120 rpm Neo/Ped: 0 rpm to 150 rpm
<b>Resolution</b>	1 rpm
<b>Accuracy</b>	Adult: 6 rpm to 120 rpm: ±2 rpm 0 rpm to 5 rpm: not specified Neo/Ped: 6 rpm to 150 rpm: ±2 rpm 0 rpm to 5 rpm: not specified
<b>Gain Selection</b>	×0.25, ×0.5, ×1, ×2, ×3, ×4, ×5
<b>Sweep</b>	6.25 mm/s, 12.5 mm/s, 25.0 mm/s, 50.0 mm/s
<b>Apnea Alarm Time Setup</b>	10 s, 15 s, 20 s (Default), 25 s, 30 s, 35 s, 40 s;

<b>NIBP</b>	
<b>Technique</b>	Oscillometry
<b>Mode</b>	Manual, Auto, Continuous, Sequence
<b>Measuring Interval in Auto Mode</b>	1/2/2.5/3/4/5/10/15/30/60/90/120/180/240/360/480 and User Define
<b>Continuous</b>	5 min, interval is 5 s
<b>Measuring Type</b>	SYS, DIA, MAP, PR
<b>Measuring Range</b>	
Adult Mode	SYS: 25 mmHg to 290 mmHg DIA: 10 mmHg to 250 mmHg MAP: 15 mmHg to 260 mmHg
Pediatric Mode	SYS: 25 mmHg to 240 mmHg DIA: 10 mmHg to 200 mmHg MAP: 15 mmHg to 215 mmHg
Neonatal Mode	SYS: 25 mmHg to 140 mmHg DIA: 10 mmHg to 115 mmHg MAP: 15 mmHg to 125 mmHg
<b>Cuff Pressure</b>	
<b>Measuring Range</b>	0 mmHg to 300 mmHg
<b>Pressure Resolution</b>	1 mmHg
<b>Maximum Mean Error</b>	±5 mmHg
<b>Maximum Standard Deviation</b>	8 mmHg
<b>Maximum Measuring Period</b>	Adult/Pediatric: 120 s Neonate: 90 s
<b>Typical Measuring Period</b>	20 s to 35 s (depend on HR/motion disturbance)
<b>Dual Independent Channel Overpressure Protection</b>	Adult: (297±3) mmHg Pediatric: (245±3) mmHg Neonatal: (147±3) mmHg

<b>EDAN Module SpO<sub>2</sub></b>	
<b>Measuring Range</b>	0% to 100%
<b>Resolution</b>	1%
<b>Data update period</b>	1 s
<b>Accuracy</b>	Adult/Pediatric: ±2% (70% to 100% SpO <sub>2</sub> ) Undefined (0% to 69% SpO <sub>2</sub> ) Neonate: ±3% (70% to 100% SpO <sub>2</sub> ) Undefined (0% to 69% SpO <sub>2</sub> )

<b>PI (Perfusion Index)</b>	
<b>Measuring Range</b>	0.00-20%, invalid PI value is -?-. 1% (10% to 20%)
<b>Resolution</b>	0.1% (1.0% to 9.9%) 0.01% (0.00% to 0.99%)

<b>Nellcor Module SpO<sub>2</sub></b>	
<b>Measuring Range</b>	1% to 100%

<b>Resolution</b>	1%
<b>Data Update Period</b>	1 s
<b>Accuracy</b>	
DS-100A, OXI-A/N (Adult)	
D-YS (Adult and Pediatric)	
OXI-P/I (Pediatric)	±3% (70% to 100% SpO <sub>2</sub> )
MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric)	±2% (70%~100% SpO <sub>2</sub> )
MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric)	±3% (60%~80% SpO <sub>2</sub> )

## PR

### PR (SpO<sub>2</sub>)

<b>Measuring range</b>	EDAN: 25 bpm to 300 bpm Nellcor: 20 bpm to 300 bpm
<b>Accuracy</b>	EDAN: ±2 bpm Nellcor: ±3 bpm (20 bpm to 250 bpm)
<b>Resolution</b>	EDAN: 1 bpm Nellcor: 1 bpm

### PR (NIBP)

<b>Measuring range</b>	EDAN: 40 bpm to 240 bpm
<b>Accuracy</b>	EDAN: ±3 bpm or 3.5%, whichever is greater
<b>Resolution</b>	EDAN: 1 bpm

### PR (IBP)

<b>Measuring range</b>	EDAN: 20 bpm to 300 bpm
<b>Accuracy</b>	EDAN: 30 bpm to 300 bpm: ±2 bpm or ±2%, whichever is greater; 20 bpm to 29 bpm: undefined
<b>Resolution</b>	EDAN: 1 bpm

## TEMP

<b>Channel</b>	2
<b>Sensor Type</b>	YSI-10K and YSI-2.252K
<b>Technique</b>	Thermal resistance
<b>Measure Parameter</b>	T1, T2, TD (the absolute value of T2 minus T1)
<b>Position</b>	Skin, oral cavity, rectum
<b>Unit</b>	°C, °F
<b>Measuring Range</b>	0°C to 50°C (32 °F to 122 °F)
<b>Resolution</b>	0.1°C (0.1 °F)
<b>Accuracy</b>	Accuracy (not including sensor): ±0.1°C Sensor accuracy: ≤ ±0.2°C
<b>Transient Response Time</b>	≤30 s

### Quick TEMP<sup>1)</sup>

<b>Measuring Range</b>	25°C ~ 45°C (monitoring mode) 35.5°C ~ 42°C (prediction mode)
<b>Sensor Type</b>	Oral/Axillary sensor, Rectal sensor
<b>Resolution</b>	0.1°C
<b>Accuracy</b>	
(not including sensor)	±0.1°C (25°C ~ 45°C) (monitoring mode)
<b>Sensor accuracy</b>	≤ ±0.2°C
<b>Update time</b>	1 s ~ 2 s
<b>Warm-up time</b>	Less than 10 seconds
<b>Prediction time</b>	Less than 30 seconds

## IBP

<b>Channel</b>	2
<b>Technique</b>	Direct invasive measurement
<b>Measuring Range</b>	Art: (0 to +300) mmHg PA/PAWP: (-6 to +120) mmHg CVP/RAP/LAP/ICP: (-10 to +40) mmHg

<b>Resolution</b>	P1/P2: (-50 to +300) mmHg
<b>Accuracy</b>	1 mmHg
<b>(not including sensor)</b>	±2% or ±1 mmHg, whichever is greater ICP: 0 mmHg to 40 mmHg: ±2 % or ±1 mmHg, whichever is greater; -10 mmHg to -1 mmHg: undefined
<b>Pressure Unit</b>	kPa, mmHg, cmH <sub>2</sub> O

### EDAN G2 Sidestream Module CO<sub>2</sub>

<b>Intended patient</b>	Adult, pediatric, neonatal
<b>Measure Parameters</b>	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR
<b>Unit</b>	mmHg, %, kPa
<b>Measuring Range</b>	EtCO <sub>2</sub> : 0 mmHg to 150 mmHg (0 % to 20%) FiCO <sub>2</sub> : 0 mmHg to 50 mmHg AwRR: 2 rpm to 150 rpm
<b>Resolution</b>	EtCO <sub>2</sub> : 1 mmHg FiCO <sub>2</sub> : 1 mmHg AwRR: 1 rpm

### EtCO<sub>2</sub> Accuracy

Typical conditions :	±2 mmHg, 0 to 40 mmHg
Ambient temperature: (25± 3) °C	±5% of reading, 41 to 70 mmHg
Barometric pressure: (760±10) mmHg	±8% of reading, 71 to 100 mmHg
Balance gas: N <sub>2</sub>	±10% of reading, 101 to 150 mmHg
Sample gas flowrate: 100 ml/min	
All conditions	±12% of reading or ±4 mmHg, whichever is greater
<b>AwRR Accuracy</b>	±1 rpm
<b>Sample Gas Flowrate</b>	70 ml/min or 100 ml/min, accuracy: ±15 ml/min
<b>Warm-up Time</b>	Display waveform within 20 s Reach the design accuracy within 2 minutes.

### Response Time

< 4 s (with 2 m gas sampling tube, sample gas flowrate: 100 ml/min)	
< 4 s (with 2 m gas sampling tube, sample gas flowrate: 70 ml/min)	

### Barometric Pressure Compensation

<b>Zero Calibration</b>	Automatic
<b>Calibration</b>	Support
<b>Apnea Alarm Delay</b>	Support
	10 s, 15 s, 20 s (Default), 25 s, 30 s, 35 s, 40 s

### Respironics Sidestream and Mainstream Modules CO<sub>2</sub>

<b>Applicable Patient Type</b>	Adult, pediatric and neonatal patients
<b>Technique</b>	Infra-red Absorption Technique
<b>Measure Parameters</b>	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR
<b>Unit</b>	mmHg, %, kpa
<b>Measuring Range</b>	EtCO <sub>2</sub> : 0 mmHg to 150 mmHg FiCO <sub>2</sub> : 3 mmHg to 50 mmHg
<b>Measuring Range</b>	AwRR: 0 rpm to 150 rpm (Mainstream) 2 rpm to 150 rpm (Sidestream)
<b>Resolution</b>	EtCO <sub>2</sub> : 1 mmHg FiCO <sub>2</sub> : 1 mmHg AwRR: 1 rpm
<b>EtCO<sub>2</sub> Accuracy</b>	±2 mmHg, 0 mmHg to 40 mmHg ±5% of reading, 41 mmHg to 70 mmHg ±8% of reading, 71 mmHg to 100 mmHg ±10% of reading, 101 mmHg to 150 mmHg ±12% of reading, RR is over 80 rpm (sidestream) There will be no degradation in performance due to respiration rate. (mainstream)
<b>AwRR</b>	±1 rpm

<b>Zero Calibration</b>	Support
<b>Apnea Alarm Delay</b>	10 s, 15 s, 20 s (Default), 25 s, 30 s, 35 s, 40 s;
<b>Barometric Pressure Compensation</b>	User setup
<b>CO<sub>2</sub> Rise Time (Mainstream)</b>	< 60 ms
<b>Sensor Response time (Sidestream)</b>	<3 seconds - includes transport time and rise time

#### Masimo Sidestream Module CO<sub>2</sub>

<b>Ambient CO<sub>2</sub></b>	≤ 800 ppm (0.08 vol%)
<b>Sampling Flow Rate</b>	(50 ± 10) sml/min
<b>Respiration Rate</b>	0 to 150 ± 1 breaths/min.
<b>Calibration</b>	No span calibration is required.
<b>Warm-up Time</b>	< 10 seconds
<b>CO<sub>2</sub> Rise Time At 50sml/min Sample Flow</b>	≤ 200 ms
<b>NomoLine ISA CO<sub>2</sub></b>	
<b>System Response Time</b>	< 3 seconds
<b>AwRR Range</b>	0 rpm to 150 rpm
<b>AwRR Accuracy</b>	± 1 rpm
<b>CO<sub>2</sub> Accuracy</b>	
Standard Conditions	(0 to 15) vol% ±(0.2 vol% + 2% of reading)
All Conditions	(15 to 25) vol% Unspecified

#### Masimo Mainstream Module CO<sub>2</sub>

<b>Respiration Rate</b>	0 to 150 ± 1 bpm.
<b>Calibration</b>	No span calibration required for the IR bench.
<b>Warm-up Time</b>	< 10 seconds (full accuracy)
<b>Rise Time (@ 10 l/min)</b>	≤ 90 ms
<b>Total System Response Time</b>	< 1 second
<b>AwRR Range</b>	0 rpm to 150 rpm
<b>AwRR Accuracy</b>	± 1 rpm
<b>CO<sub>2</sub> Accuracy</b>	
Standard Conditions	(0 to 15) vol% ±(0.2 vol% + 2% of reading)
All Conditions	±(0.3 kPa + 4% of reading)

#### Safety Specifications

<b>Compliant with Standards</b>	IEC 60601-1: 2005+A1 :2012; IEC 60601-1-2: 2014; EN 60601-1: 2006+A1 :2013; EN 60601-1-2: 2015; IEC 80601-2-49: 2018
<b>Anti-electroshock type</b>	Class I equipment and internal powered equipment
<b>Anti-electroshock degree</b>	CF: ECG (RESP), TEMP, IBP, Quick Temp BF: SpO <sub>2</sub> , NIBP, CO <sub>2</sub>
<b>Ingress Protection</b>	IPX1

#### Environmental Specifications

<b>Temperature</b>	Working : +0°C to +40°C (32°F ~ 104°F) Transport and Storage: -20°C to +55°C (-4°F ~ 131°F)
<b>Humidity</b>	Working: 15%RH to 95%RH (non-condensing) Transport and Storage: 15%RH to 95%RH (non-condensing)
<b>Altitude</b>	Working: 86 kPa to 106 kPa Transport and Storage: 70 kPa to 106 kPa